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
ROLLSTÄNDER LIGHT DUTY

ASSEMBLY INSTRUCTION

CART LIGHT DUTY

Name: CIM Rollständer Light Duty

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PLEASE READ THIS INSTRUCTION CAREFULLY BEFORE INSTALLING!

Clarification of non-standardized Symbols:

-  Warning of a danger point
-  Important information
-  Danger of squeezing fingers

AREA OF APPLICATION/ INTENDED USE

CIM medical products are intended to mount and position medical electrical equipment (hereafter referred to as ME equipment) and accessories to those and are used in operating theatres as well as treatment rooms as auxiliary devices for doctors and medical staff.

GENERAL NOTES

Verify that your cart configuration is approved by the CIM med GmbH.

WARNING Installation, Operation and Maintenance 

This section contains warnings regarding the installation, operation and maintenance of the cart. This section must be read in its entirety before installing, operating and maintaining the cart. Failure to follow these warnings may result in damage to equipment or injury to people.

Installation, operation and maintenance may only be carried out by qualified personnel. Cable integration and assembly may only be done by qualified technical personnel. CIM med GmbH cannot be held responsible for any damage due to poor installation.

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GENERAL NOTES

Please consider carefully the ME equipment being mounted. CIM med GmbH recommends that the hospital's risk management personnel verify that the application is appropriate prior to installation and use of the cart.

Please verify the weight added. The maximum load, marked on a cart or in the accompanying drawing must not be exceeded. Too high weight load can lead to tilting of a cart at surfaces with 10° incline. If you ordered a cart with counterweight make sure that there is one at the base of your cart. Check that the counterweight is properly installed before you start loading.

An adjustment of the payload requires an approval by CIM med GmbH engineering design department.

Product modifications, disassembling and reassembling and/or mounting any of the components in any other way than described below are permitted only after written approval by the CIM med GmbH construction and engineering department. CIM med GmbH cannot be held responsible for any damage or injury caused by product modifications and/or amendments without prior permission by CIM med GmbH.

⚠️ WARNING: MORTAL DANGER BY ELECTRIC SHOCK

The electric cables of ME equipment and ME systems may only be connected to the power source and to the medical device respectively after finishing installation. Prior to maintenance or disassembly all power cables must be disconnected from the power source. Only cables with double insulation may be used for running through the support arm.

TOOLS REQUIRED

Torx wrench TX 45 for column
Philips screwdriver PH2 for counterweight
Open-end wrench size 21 for castors

CLEANING

CIM carts may be cleaned with mild, non-abrasive solutions commonly used in hospital environment. Do not use steel wool or other abrasive materials. Avoid all liquids entering the cart. Your cart has been designed to the highest specifications to make cleaning as easy as possible and to avoid liquids from penetrating. CIM med GmbH makes no claims concerning the efficacy of the disinfection processes. It is recommended to test any cleaning solution on a small area of the cart not visible to verify compatibility. Please consult your hospital's infection control manager.

MAINTENANCE: MANNER AND FREQUENCY

Very little maintenance is required during the life span of the cart. Please check all screws on a regular basis (at least once a year) and tighten them if necessary.

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In case of damage done to the load-bearing parts of the system it has to be strictly refrained from further usage of the cart.

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Inspection must include the following steps:

Grasp the ME equipment and tilt it forward and back, through its entire range of motion. There should be enough tension or resistance in the tilt mechanism to prevent the ME equipment from tilting forward unexpectedly. If necessary, the tilt tension may be adjusted in accordance with section „Tilt adjustment“.

Check that the castors are firmly bolted to the base. Move the cart in a circle, making sure that all castors rotate smoothly and turn around the axis. Check the functionality of the brakes.

Failure to periodically inspect and adjust the cart as instructed may result in damage to equipment or injury.

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SAFEKEEPING

For any questions that may arise at a later time please safekeep our manual!

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
MONITOR ADAPTATIONS

REAR-MOUNTED DEVICES

For rear-mounted devices, see separate assembly instructions:

VESA 75/ 100

- ➔ Assembly Instruction Monitor Adaptation VESA 75, VESA 75/ 100

VESA 75

- ➔ Assembly Instruction Monitor Adaptation VESA 75, VESA 75/ 100

VESA 75 QUICK RELEASE

- ➔ Assembly Instruction Monitor Adaptation Quick Release VESA 75, VESA 75/ 100

VESA 75/100 QUICK RELEASE

- ➔ Assembly Instruction Monitor Adaptation Quick Release VESA 75, VESA 75/ 100

BOTTOM-MOUNTED DEVICES

For bottom-mounted devices, see separate assembly instructions:

5" PLUNGE PLATE

- ➔ Assembly Instruction Monitor Adaptation 5" plunge plate

MINDRAY QUICK RELEASE

- ➔ Assembly Instruction Monitor Adaptation Mindray Quick Release

PHILIPS TABLE TOP MOUNT

- ➔ Assembly Instruction Monitor Adaptation Philips Table Top Mount

SPACELABS XPREZZON BACKPACK

- ➔ Assembly Instruction Monitor Adaptation Spacelabs XPREZZON Backpack

FUKUDA MOUNTING ADAPTER

- ➔ Assembly Instruction Monitor Adaptation Fukuda

NOTEBOOK TRAY


- ➔ Assembly Instruction Notebook Tray for height adjustable arm

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

ASSEMBLY STEPS


ASSEMBLY STEPS

1. Loosen the (5) M5x12 pan head screws and remove the counterweight if existing.

ⓘ NOTE: Hold the counterweight while loosening the screws.
2. Push the accessory on the designated side into the column and mount it to the given height. Please refer to the corresponding assembly instruction of the accessory.


3. Mount the column to the base with (4) M8x35 cylinder head bolts and (4) contact washers. Tighten the screws securely.
4. Mount the counterweight by using the (5) M5x12 pan head screws.

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


ASSEMBLY STEPS

5. Choose the required length of the plastic covers and then cut them.

6. Clip the plastic covers into the column.

TILT ADJUSTMENT

7. To adjust the overall tilt tension, tighten or loosen the both hex screws with the hex key 4 mm provided. By tightening both of the screws the tilt function can be locked completely. Maximum permissible tightening torque is 8 Nm (70 in-lbs).

⚠ WARNING: In case of overtightening the screw (> 10 Nm/ 88 in-lbs) the casting respectively the threads may break which leads to loss of function.

⚠ WARNING: When loosening both hex screws the ME equipment will tilt forward.

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MAINTENANCE PLAN

CIM med GmbH recommends the following maintenance schedule for annual inspection.

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Die CIM med GmbH empfiehlt folgenden Wartungsplan zur jährlichen Kontrolle.