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

MONITORANBINDUNG AN ONDAL ACROBAT 2000/ 3000/ 5000

ASSEMBLY INSTRUCTION

MONITOR ADAPTATION AT ONDAL ACROBAT 2000/ 3000/ 5000

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PLEASE READ THIS INSTRUCTION CAREFULLY BEFORE INSTALLING!

Clarification of non-standardized Symbols:

 Warning of a danger point

 Important information

AREA OF APPLICATION/ INTENDED USE

Cable Integrated Mounting systems (CIM) are intended to mount and position monitors, mouse/keyboard, laptops and other medical devices or equipment (hereafter referred to as ME equipment). They are mounted to walls, carts, tables, other surfaces or machines using specific mounting adapters. CIM systems do not fulfil any specific medical purpose on their own.

GENERAL NOTES

Verify that the Ondal arm has been installed and approved according to installation instructions. Please verify that the total weight do not exceed the permissible total weight of your Ondal arm. The permissible total weight is mentioned on the corresponding data sheet.

WARNING Installation, Operation and Maintenance 

This section contains warnings regarding the installation, operation and maintenance of the monitor adaptation at Ondal Acrobat 2000/ 3000/ 5000 (hereafter referred to as monitor adaptation). This section must be read in its entirety before installing, operating and maintaining the monitor adaptation. Failure to follow these warnings may result in damage to equipment or injury to people.

Installation, operation and maintenance may only be carried out by qualified personnel. Cable integration and assembly may only be done by qualified technical personnel. CIM med GmbH cannot be held responsible for any damage due to poor installation.

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GENERAL NOTES

Please consider carefully the ME equipment being mounted and the proximity of the mounting assembly to other equipment, hospital personnel, and the patient. CIM med GmbH recommends that the hospital's risk management personnel verify that the application is appropriate prior to installation and use of the monitor adaptation.

Please check the weight to be mounted. The weight of the monitor adaptation and the mounted ME equipment must not exceed the permissible total weight of the Ondal arm.

Product modifications, disassembling and reassembling and/or mounting any of the components in any other way than described below are permitted only after written approval by the CIM med GmbH construction and engineering department. CIM med GmbH cannot be held responsible for any damage or injury caused by product modifications and/ or amendments without prior permission by CIM med GmbH.

⚠️ WARNING: MORTAL DANGER BY ELECTRIC SHOCK

The electric cables of ME equipment and ME systems may only be connected to the power source and to the medical device respectively after finishing installation. Prior to maintenance or disassembly all power cables must be disconnected from the power source. Only cables with double insulation may be used.

⚠️ WARNING: EARTHING CONNECTION

The monitor adaptation is delivered with components' grounding. It is the healthcare institution's responsibility to connect the Ondal fixing to the electrical grounding of the Ondal arm.

TOOLS REQUIRED

hex key 2.5 for alu cover cap

CLEANING AND DISINFECTION

The monitor adaptation may be cleaned with mild, non-abrasive solutions commonly used in hospital environment. Do not use steel wool or other abrasive materials. Avoid all liquids entering the monitor adaptation. The monitor adaptation has been designed to the highest specifications to make cleaning as easy as possible and to avoid liquids from penetrating. CIM med GmbH makes no claims concerning the efficacy of the disinfection processes. It is recommended to test any cleaning solution on a small area not visible to verify compatibility. Please consult your hospital's infection control manager.

Exemplary disinfection tests have been carried out with the following active substance groups: Alcoholic quick disinfectants (e.g. Bacillol AF/Hartmann), aldehyde and quaternary ammonium compounds (e.g. B5/Orochemie), alcohol-free quaternary ammonium compounds (e.g. Sani-Cloth AF3/PDI Healthcare). These compounds are VAH / DGHM listed and EPA registered and must be used as per the manufacturers' instructions. The surface finish will be permanently damaged by acetone and trichloroethylene.

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GENERAL NOTES

MAINTENANCE: MANNER AND FREQUENCY

Very little maintenance is required during the life span of the Ondal fixing; however, it is recommended to check all pivot/ swivel points as well as the bearing point of the Ondal support arm regularly. Please check all screws on a regular basis (at least once a year) and tighten them if necessary. In case of damage done to the load-bearing parts of the system it has to be strictly refrained from further usage of the monitor adaptation.

Inspection must include the following steps:

Grasp the ME equipment and tilt it forward and back, through its entire range of motion. There should be enough tension or resistance in the tilt mechanism to prevent the ME equipment from tilting forward unexpectedly. If necessary, the tilt tension may be adjusted in accordance with section „Tilt adjustment“.

Rotate the Ondal fixing around its axis. There should be no resistance in the range of motion. Check the functionality of its limit of rotation.

REPORTING OBLIGATION

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SAFEKEEPING

For any questions that may arise at a later time please safekeep our manual!

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

ASSEMBLY STEPS


1. Loosen the (4) M4x8 pan head screws and remove the alu cover cap.

2. Feed the cables one after the other through the column.


NOTE: Always start with the largest cable first if you have multiple cables.

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
3. Feed cables into the horizontal column from below before feeding them through the column to the monitor adaptation.

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

4. Exit the cables below the swivel part.

5. Connect the components' grounding of the monitor mount with the Ondal arm by using the shown M4x12 pan head screw.

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


ASSEMBLY STEPS

6. Mount the alu cover cap using the (4) M4x8 pan head screws.
7. Define the necessary length of the channel cover and cut it.

8. Clip the channel cover into the column.
9. Mount the Ondal fixing to the Ondal arm. Please refer to the corresponding assembly instruction of the Ondal arm.

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CIM med GmbH recommends the following maintenance schedule for annual inspection.

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Die CIM med GmbH empfiehlt folgenden Wartungsplan zur jährlichen Kontrolle.