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
Clarification of non-standardized Symbols:

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**AREA OF APPLICATION/ INTENDED USE**

Cable Integrated Mounting systems (CIM) are intended to mount and position monitors, mouse/keyboard, laptops and other medical devices or equipment (hereafter referred to as ME equipment). They are mounted to walls, carts, tables, other surfaces or machines using specific mounting adapters. CIM systems do not fulfil any specific medical purpose on their own.

**GENERAL NOTES**

Make sure that the channel into which you are inserting single mount (hereinafter referred to as the support arm), is correctly mounted and can withstand the weight load. Please check the weight of the ME equipment to be mounted. The permissible total weight is mentioned on the corresponding data sheet.

**WARNING Installation, Operation and Maintenance** 

This section contains warnings regarding the installation, operation and maintenance of the support arm. This section must be read in its entirety before installing, operating and maintaining the support arm. Failure to follow these warnings may result in damage to equipment or injury to people.

Installation, operation and maintenance may only be carried out by qualified personnel. Cable integration and assembly may only be done by qualified technical personnel. CIM med GmbH cannot be held responsible for any damage due to poor installation.

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## GENERAL NOTES

Do not position the support arm above a patient. Note that the support arm has a wide range of motion both up/down and side to side. Please consider carefully the ME equipment being mounted and the proximity of the mounting assembly to other equipment, hospital personnel, and the patient. CIM med GmbH recommends that the hospital's risk management personnel verify that the application is appropriate prior to installation and use of the support arm.

Ensure that the weight of the instrument being mounted does not exceed the maximum load indication.

Product modifications, disassembling and reassembling and/or mounting any of the components in any other way than described below are permitted only after written approval by the CIM med GmbH construction and engineering department. CIM med GmbH cannot be held responsible for any damage or injury caused by product modifications and/or amendments without prior permission by CIM med GmbH.

### **⚠️ WARNING: MORTAL DANGER BY ELECTRIC SHOCK**

The electric cables of ME equipment and ME systems may only be connected to the power source and to the medical device respectively after finishing installation. Prior to maintenance or disassembly all power cables must be disconnected from the power source. Only cables with double insulation may be used for running through the support arm.

### **⚠️ WARNING: EARTHING CONNECTION**

The support arm is delivered with components' grounding. It is the health-care institution's responsibility to control the grounding of the rails. The support arm is connected to the channel via the grounding screw.

### **TOOLS REQUIRED**

Torx wrench TX 30 for sliding block, support plate  
 Torx wrench TX 45 for support arm  
 Hex key 3.0 for grounding screw  
 Hex key 5.0 for sliding block  
 Hex key 6.0 for pole clamp

### **CLEANING AND DISINFECTION**

CIM support arms may be cleaned with mild, non-abrasive solutions commonly used in hospital environment. Do not use steel wool or other abrasive materials. Avoid all liquids entering the support arm. Your support arm has been designed to the highest specifications to make cleaning as easy as possible and to avoid liquids from penetrating. CIM med GmbH makes no claims concerning the efficacy of the disinfection processes. It is recommended to test any cleaning solution on a small area of the support arm not visible to verify compatibility. Please consult your hospital's infection control manager.

Exemplary disinfection tests have been carried out with the following active substance groups: alcoholic quick disinfectants (e.g. Bacillol AF/Hartmann), aldehyde and quaternary ammonium compounds (e.g. B5/Orochemie), alcohol-free quaternary ammonium compounds (e.g. Sani-Cloth AF3/PDI

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CIM med GmbH  
Margot-Kalinke-Str. 9, 80939 München

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## GENERAL NOTES

Healthcare). These compounds are VAH / DGHM listed and EPA registered and must be used as per the manufacturers' instructions. The surface finish will be permanently damaged by acetone and trichloroethylene.

### MAINTENANCE: MANNER AND FREQUENCY

Very little maintenance is required during the life span of the support arm; however, it is recommended to check all pivot/ swivel points as well as all bearing points regularly. Please check all screws on a regular basis (at least once a year) and tighten them if necessary.

The support arm is designed for 26,000 movements or 13,000 cycles (one cycle = a full movement from one extreme to the other and back again). A daily number of 10 movements on 365 days (= operation on 7 days a week) means 3,650 movements per year which correspond to a durability of 7.1 years.

In case of damage done to the load-bearing parts of the system it has to be strictly refrained from further usage of the mounting system.

### Inspection must include the following steps:

Grasp the support arm and pivot it from side to side at the attachment. The support arm should pivot with some tension or resistance.

Failure to periodically inspect and adjust the support arm as instructed may result in damage to equipment or injury.

### REPORTING OBLIGATION

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### SAFEKEEPING

For any questions that may arise at a later time please safekeep our manual!

### ASSISTANCE REQUIRED

If you are unable to adjust the arm or require service, contact CIM med technical support: Phone +49 89 978 94 080 or your local dealer.

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


# ASSEMBLY STEPS

1. Loosen the (2) M6x25 pan head screws for one turn.
2. Insert (2) cylinder pins into the lower holes. Mount support arm with (2) M8x30 cylinder head screws and (2) contact washers.  
Torque: 25 Nm


**⚠ WARNING:** The use of contact washers is absolutely necessary for a safe grounding of components.

3. Thread single mount into the channel of anaesthesia machine.

**ⓘ NOTE:** Loosen the screw connection of the T-nuts to allow threading.


4. Tighten the (2) M6x20 cylinder head screws and (2) M6x25 pan head screws.

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

## ASSEMBLY STEPS


5. Fasten the (2) M6x25 pan head screws to fix the bridge and support plate.
6. Mount the (1) M6x20 threaded pin.

**⚠ WARNING:** The use of threaded pin is absolutely necessary for a safe grounding of components.

7. Remove the (3) cover caps and loosen the (2) M8x55 cylinder head screws. Insert the infusion pole.
8. First tighten the (2) M8x55 cylinder head screws evenly. Now tighten the (1) M6x8 threaded pin and cover all screws with the (3) cover caps.

**⚠ WARNING:** The use of threaded pin is absolutely necessary for a safe grounding of components.

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CIM med GmbH recommends the following maintenance schedule for annual inspection.

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Die CIM med GmbH empfiehlt folgenden Wartungsplan zur jährlichen Kontrolle.