


Cable Management - Cable hook



Illustration	Description	Fixing	Part no.
	Down post with cable hook To mount to front joint of arm	To mount to non-height adjustable arms Length 150 mm / 5.9" Length 300 mm / 11.8"	Z-KH6-70.1_150 Z-KH6-70.1_300
		To mount to height adjustable arms Length 150 mm / 5.9" Length 300 mm / 11.8"	Z-KH6-71.1_150 Z-KH6-71.1_300

General notes

Custom-designed solutions and MRI applications

Please contact us for **custom-designed** solutions.

Product marking

All products which are classified as medical product risk class 1 are CE marked. As **medical product risk class 1** they are conform with the regulations (EU) 2017/745 (MDR) and IEC 60601-1. Other standards of the IEC 60601 series are not applicable to our products, regardless of their version.

MD A stand-alone product that is already classified as a medical product class 1 before being integrated into an overall medical device. The product is CE-marked. A declaration of conformity is available.

~~**MD**~~ Accessories or spare parts that have no independent function and are therefore not subject to the MDR. These products are not CE-marked. A declaration of conformity is not required.

Warranty

All CIM arms have a 5-year warranty to be free of defects in material and functionality from the date of delivery.

