


## Wall channel

Illustration	Description	Length	Part no.
	<p><b>CIM wall channel</b></p> <p>Vertical wall channel including channel cover and end caps</p> <p>Individual lengths on request!</p> <p>Channel cover available in two colours (C): RAL 7024 graphite grey (gr) or RAL 9016 traffic white (w)</p>	<p>480 mm / 18.9"</p> <p>720 mm / 28.3"</p> <p>960 mm / 37.8"</p> <p>1200 mm / 47.2"</p>	<p>Z-WK.C_480</p> <p>Z-WK.C_720</p> <p>Z-WK.C_960</p> <p>Z-WK.C_1200</p>

## General notes

---

### Product marking

All products which are classified as medical product risk class 1 are CE marked. As **medical product risk class 1** they are conform with the regulations (EU) 2017/745 (MDR) and IEC 60601-1. Other standards of the IEC 60601 series are not applicable to our products, regardless of their version.

### Warranty

All CIM arms have a 5-year warranty to be free of defects in material and functionality from the date of delivery.

### Custom-designed solutions and MRI applications

Please contact us for **custom-designed** solutions.

