

## ASSEMBLY GADGET FOR CABLING A CIM ARM

In order to facilitate the cabling of the support arms on site and, above all, to save assembly time in the hospital, we have designed an assembly gadget that enables pre-assembly in hospital and operating rooms. This means that all support arms with internal cables can be pre-assembled outside the hospital rooms / operating rooms and valuable time in the final assembly is saved.

Z-AG-01



### MOUNTING CASE INCLUDES

- 1x Transport case
- 2x Bar clamps
- 1x Base plate
- 1x Pole
- 1x Column
- 1x Adapter horizontal rail
- 1x Adapter vertical rail

### FIXINGS

- 1V Universal pole clamp for pole  $\varnothing$  23 - 40 mm / 0.9 - 1.6"
- 1C Fixed pole clamp for pole  $\varnothing$  38 mm / 1.5"
- 2F Horizontal rail clamp for European rail 10 x 25 mm
- 3 Vertical rail clamp for European rail 10 x 25 mm
- 4S Adapter for vertical wall channel

