


# Mindray monitor-specific mounting adapter



Type of monitor	Description	Part No.
		<del>MD</del>
PM 8000	Mounting adapter for 5" plunge plate	MD-AP-PM8
BeneView T8 / BeneView T9	Mounting adapter for 5" plunge plate	MD-AP-T8
BeneView T5 / PM 9000	Mounting adapter for 5" plunge plate	MD-AP-T5
BeneVision N15 / N17	Mounting adapter for 5" plunge plate	MD-AP-N
<p>Adapter plates are only supplied with CIM support arms, counters tops or carts. Exceptions can be made on request, however, please note that price and part number will change when sold separately (K-MD-AP-...).</p>		