


Handle



Illustration	Description	Notes	Weight	Part no.
	Dual handle with VESA fixing	For monitors with a maximum width of 430 mm / 16.9"	1.5 kg / 3.3 lbs	Z-HDG.1_43
		570 mm / 22.4"	1.9 kg / 4.2 lbs	Z-HDG.1_57

Example of use Z-HDG.1



General notes

Custom-designed solutions and MRI applications

Please contact us for **custom-designed** solutions.

Product marking

All products which are classified as medical product risk class 1 are CE marked. As **medical product risk class 1** they are conform with the regulations (EU) 2017/745 (MDR) and IEC 60601-1. Other standards of the IEC 60601 series are not applicable to our products, regardless of their version.

MD A stand-alone product that is already classified as a medical product class 1 before being integrated into an overall medical device. The product is CE-marked. A declaration of conformity is available.

~~MD~~ Accessories or spare parts that have no independent function and are therefore not subject to the MDR. These products are not CE-marked. A declaration of conformity is not required.

Warranty

All CIM products have a 5-year warranty to be free of defects in material and functionality from the date of delivery.

How to find the right product for your needs?

No problem with our configurator

<https://cim-configurator.com>

