


## Cable Management - Cable hook



Illustration	Description	Fixing	Weight	Part no.	
	<b>Down post with cable hook</b> To mount to front joint of arm	To mount to non-height adjustable arms Length 150 mm / 5.9"	0.3 kg / 0.7 lbs	Z-KH6-70.1_150	
		Length 300 mm / 11.8"	0.4 kg / 0.9 lbs	Z-KH6-70.1_300	
			To mount to height adjustable arms Length 150 mm / 5.9"	0.3 kg / 0.7 lbs	Z-KH6-71.1_150
			Length 300 mm / 11.8"	0.4 kg / 0.9 lbs	Z-KH6-71.1_300

## General notes

---

### Custom-designed solutions and MRI applications

Please contact us for **custom-designed** solutions.

### Product marking

All products which are classified as medical product risk class 1 are CE marked. As **medical product risk class 1** they are conform with the regulations (EU) 2017/745 (MDR) and IEC 60601-1. Other standards of the IEC 60601 series are not applicable to our products, regardless of their version.

**MD** A stand-alone product that is already classified as a medical product class 1 before being integrated into an overall medical device. The product is CE-marked. A declaration of conformity is available.

~~**MD**~~ Accessories or spare parts that have no independent function and are therefore not subject to the MDR. These products are not CE-marked. A declaration of conformity is not required.

### Warranty

All CIM products have a 5-year warranty to be free of defects in material and functionality from the date of delivery.

---

## How to find the right product for your needs?

No problem with our configurator

<https://cim-configurator.com>

