


# European rail



Illustration	Description	Length	Weight	Part no.
	<p><b>European rail 10 x 25 mm</b></p> <p>To mount to CIM wall channel including adapter for vertical wall channel</p>	<p>300 mm / 11.8"</p> <p>400 mm / 15.7"</p> <p>500 mm / 19.7"</p>	<p>0.4 kg / 0.9 lbs</p> <p>0.5 kg / 1.1 lbs</p> <p>0.5 kg / 1.1 lbs</p>	<p>Z-NS-4Z.1_300</p> <p>Z-NS-4Z.1_400</p> <p>Z-NS-4Z.1_500</p>

## General notes

---

### Custom-designed solutions and MRI applications

Please contact us for **custom-designed** solutions.

### Product marking

All products which are classified as medical product risk class 1 are CE marked. As **medical product risk class 1** they are conform with the regulations (EU) 2017/745 (MDR) and IEC 60601-1. Other standards of the IEC 60601 series are not applicable to our products, regardless of their version.

**MD** A stand-alone product that is already classified as a medical product class 1 before being integrated into an overall medical device. The product is CE-marked. A declaration of conformity is available.

~~MD~~ Accessories or spare parts that have no independent function and are therefore not subject to the MDR. These products are not CE-marked. A declaration of conformity is not required.

### Warranty

All CIM products have a 5-year warranty to be free of defects in material and functionality from the date of delivery.

---

## How to find the right product for your needs?

No problem with our configurator

<https://cim-configurator.com>

