

# BASIC MOBILE STAND

SIMPLE, MOBILE, RELIABLE



Compact, agile, and built to last: The Basic Mobile Stand is the perfect solution for everyday clinical use when only one device is required. With effortless handling, smooth maneuverability, and easy-to-clean surfaces, it ensures reliable performance in daily routines. Every configuration is fully tested and safety-approved, and with a 5-year warranty, it delivers lasting quality and peace of mind.

## KEY FEATURES

- + Versatile use**  
Ideal as a monitor cart or for patient monitoring.
- + Lightweight & maneuverable**  
Optimal mobility, even in confined spaces.
- + Compact design**  
1158 x 500 x 559 mm (H x W x D, without accessories)
- + Certified & safe**  
Compliant with MDR (EU 2017/745, Class I) and IEC 60601-1.
- + Strong & stable**  
Supports up to 8 kg on the monitor mount and 5 kg in the basket.


## ACCESSORIES

Delivered with a handle as standard. Additional practical accessories are available to adapt the cart perfectly to your individual requirements.



## MONITOR PRESENTATIONS

Which monitor would you like to mount?  
 CIM med offers solutions for a wide range of medical device manufacturers.



## TECHNICAL DATA

