

## Infusion pump mount - accessories



### Material

Stainless steel

### Part number

Infusion pole ø 25 mm / 1"

Part number	Length	Product weight
Z-IF_600	600 mm / 23.6"	0.7 kg / 1.5 lbs
Z-IF_1200	1200 mm / 47.2"	1.3 kg / 2.9 lbs

Infusion pole ø 25 mm / 1" with infusion bag hooks

Part number	Length	Product weight
Z-IF6_600	600 mm / 23.6"	0.9 kg / 2.0 lbs
Z-IF6_1200	1200 mm / 47.2"	1.5 kg / 3.3 lbs

## General notes

---

### Custom-designed solutions and MRI applications

Please contact us for **custom-designed** solutions.

### Product marking

All products which are classified as medical product risk class 1 are CE marked. As **medical product risk class 1** they are conform with the regulations (EU) 2017/745 (MDR) and IEC 60601-1. Other standards of the IEC 60601 series are not applicable to our products, regardless of their version.

**MD** A stand-alone product that is already classified as a medical product class 1 before being integrated into an overall medical device. The product is CE-marked. A declaration of conformity is available.

~~**MD**~~ Accessories or spare parts that have no independent function and are therefore not subject to the MDR. These products are not CE-marked. A declaration of conformity is not required.

### Warranty

All CIM products have a 5-year warranty to be free of defects in material and functionality from the date of delivery.

---

## How to find the right product for your needs?

No problem with our configurator

<https://cim-configurator.com>

