

CART ULTRA-LIGHT DUTY

SIMPLE, MOBILE, RELIABLE

Compact, agile, and built to last: The ultra-light duty cart is the perfect solution for everyday clinical use when only one device is required. With effortless handling, smooth maneuverability, and easy-to-clean surfaces, it ensures reliable performance in daily routines. Every configuration is fully tested and safety-approved, and with a 5-year warranty, it delivers lasting quality and peace of mind.

KEY FEATURES

-  **Versatile use**
Ideal as a monitor cart or for patient monitoring.
-  **Lightweight & maneuverable**
Optimal mobility, even in confined spaces.
-  **Compact design**
1158 x 500 x 559 mm (H x W x D, without accessories)
-  **Certified & safe**
Compliant with MDR (EU 2017/745, Class I) and IEC 60601-1.
-  **Strong & stable**
Supports up to 8 kg on the monitor mount and 5 kg in the basket.

PART NUMBER

R-003B.C1 **Cart Ultra-Light Duty**

ACCESSORIES

Additional practical accessories are available to adapt the cart perfectly to your individual requirements.

Dual Cable Hook



Triple Cable Hook



Basket medium



Basket small



Handle



MONITOR PRESENTATIONS

Which monitor would you like to mount?

CIM med offers solutions for a wide range of medical device manufacturers.



TECHNICAL DATA

