

Fixings for Frame F2 / F5 / Canvas

Illustration	Description	Part no.	
	L-Bracket for Frame F2 / F5 / Canvas Without fixing	GE-MB-0.1	MD
	L-Bracket for Frame F2 / F5 / Canvas Fixing: fixed pole clamp to mount to pole ø 38 mm / 1.5"	GE-MB-1C.1	MD
	L-Bracket for Frame F2 / F5 / Canvas Fixing: adapter for vertical wall channel H 80 mm / 3.2"	GE-MB-4.1	MD
	L-Bracket for Frame F2 / F5 / Canvas Fixing: horizontal rail clamp to mount to horizontal European rails 10 x 25, 8 x 35, 10 x 30, 10 x 50 mm	GE-MB-6L.1	MD
	C-Bracket for Frame F2 / F5 / Canvas L 150 mm / 5.9" For non-height adjustable arms For height adjustable arms	GE-MB-DPR-70.1_150 GE-MB-DPR-71.1_150	MD