

Accessories for CARESCAPE ONE

Illustration	Description	Notes	Part no.	
	Down post, rotatable with cable hook and GE Carescape ONE adaption	For non-height adjustable arms: L 300 mm / 11.8", max. load 10 kg / 22 lbs For height adjustable arms: L 300 mm / 11.8", max. load 10 kg / 22 lbs	Z-DPR-GE-70.C1_300 Z-DPR-GE-71.C1_300	MD
	Down post, rotatable with cable hook	For non-height adjustable arms: L 300 mm / 11.8", max. load 10 kg / 22 lbs For height adjustable arms: L 300 mm / 11.8", max. load 10 kg / 22 lbs	Z-KH6-70.1_300 Z-KH6-71.1_300	MD
	Adapter for CARESCAPE ONE	Mounting to down post	Z-GE-1.1	MD